

reGO X-ray and ultrasound accessories

reGO - Kassettenhalterung Typ R für Rasterwandstative

-cassette holder type R for grid wall tripods



Die R-Halterung eignet sich ausgezeichnet für Aufnahmen ohne Laufraster. Sie ist aus hochwertigem, eloxiertem Aluminium gefertigt und zeichnet sich durch hohe Stabilität, verbunden mit einem geringen Gewicht von ca. 2,5 kg aus.

- Befestigung durch Einhängen am Rasterwandstativ,
- passend für alle gängigen Frontplatten mit Flachrastschiene, DIN-Profilschiene,
- lieferbar für verschiedene Kassettenformate

The cassette holder is recommended for X-rays without sliding grids. It's made of high-quality anodized aluminium which guarantees a high grade of stability in spite of its low weight of only approx. 2.5 kg.

Fixed by hanging onto the grid wall tripod.

Suitable for all standard front plates with profile sliding mechanism.

Available for various cassette sizes.

Art.Nr. 102111 für Kassetten:
for cassettes : 30 x 90 cm, 20 x 60, 20 x 80, 30 x 60, 30 x 80 cm

Art.Nr. 102112 für Kassetten:
for cassettes: 30 x 120 cm, 20 x 96 cm

Art.Nr. 102113 für Kassetten:
for cassettes: 13 x 18 bis/to 35 x 43 cm