

(D)

Rego-Raster ab 60 lp/cm sind auch für das digitale Röntgen geeignet, bei der Verwendung von Kassetten mit Speicherfolien von z.B. Fuji, Agfa und andere. Sie können als Rasterbrücken oder als Außenraster eingesetzt werden.

Große Außenraster (z.B. 35x84, 35x124 cm) zum Erstellen von Aufnahmen am Vertikalstativ sind rundum mit einem Kohlefaserrahmen eingefasst und strahlenseitig von einer Kohlefaserabdeckung geschützt.

Die Raster werden vor die Kassette in die Kassettenhalterung gestellt, welche am Vertikalstativ eingehängt wird. (Abb.)

(GB)

Rego-grids from 60 lp/cm are also suitable for digital X-Ray when using cassettes with imaging plates from Fuji and Kodak, among others.

They can be used as slip-on grids or outside grids. Large outside grids (e.g. 35 x 84, 35 x 124 cm) for the taking of images on a vertical stand are installed in carbon fibre frames, the tube side being protected by a carbon fibre covering.

The grids are placed in front of the cassettes in the cassette holder, this then being placed in the vertical stand (see illustration).



Bitte betrachten Sie die Angaben und Abbildungen in dem Prospekt als unverbindlich.
Wir behalten uns Änderungen des Liefergegenstandes gegenüber den hier gezeigten Modellen vor.
All prints in this brochure are not binding and subject to alterations.

Les indications et figures dans ce prospectus s'entendent sans engagement. Nous réservons le droit d'y apporter des modifications



DIN EN ISO 9001 : 2008
REG.NR.: 292498 QM08

rego X-ray GmbH
X-RAY ACCESSORIES & MECHANICAL ENGINEERING
Stuttgarter Str. 5, D-86154 Augsburg
Tel. 0049/821-450435-0 Fax 0049/821-450435-11
e-mail: info@rego.de <http://www.rego.de>

